




Recombinant Human Tyrosine-protein kinase SYK (SYK)

Product Code	CSB-EP022998HU
Abbreviation	Recombinant Human SYK protein
Uniprot No.	P43405
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.72.
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	<pre>MASSGMADSANHLPFFFGNITREEAEDYLVQGGMSDGLYLLRQSRNYLGGFA LSVAHGRKAHHYTIERELNGTYAIAGGRTHASPADLCHYHSQESDGLVCLLKK PFNRPQGVQPKTGFEDLKENLIREYVKQTNLQGGQALEQAIISQKPQLEKLI TTAHEKMPWFHKGISREESEQIVLIGSKTNGKFLIRARDNNGSYALCCLLHEGKV LHYRIDKDKTGKLSIPEGKKFDLWQLVEHYSYKADGLLRVLTVPCQKIGTQGN VNFGRPQLPGSHPATWSAGGIISRIKSYSFPKPGHRKSSPAQGNRQUESTVS FNPYEPELAPWAADKGPQREALPMDTEVYESPYADPEEIRPKEVYLDRKLLTL EDKELGSGNFGTVKKGYYQMKKVVKTVAVKILKNEANDPALKDELLAEANVM QQLDNPYIVRMIGICEAESWMLVMEAEELGPLNKYLQQNRHVKDKNIIELVHQ VSMGMKYLEESNFVHRDLAARNVLLVTQHYAKISDFGLSKALRADENYYKAQT HGKWPVKWYAPECINYYKFSSKSDVWSFGVLMWEAFSYGQKPYRGMKGSE VTAMLEKGERMGCPAGCPREMYDLMNLCWTDVENRPGFAAVELRLRNYYY DVVN</pre>
Research Area	Immunology
Source	E.coli
Target Names	SYK
Expression Region	1-635aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	79.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.