




# Recombinant Human Osteopontin (SPP1), partial

<b>Product Code</b>	CSB-EP022603HU2
<b>Abbreviation</b>	Recombinant Human SPP1 protein, partial
<b>Uniprot No.</b>	P10451
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.81.
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	IPVKQADSGSSEEKQLYNKYPDAVATWLNPDPSQKQNLLAPQNAVSSEETND FKQETLPSKSNESHDMDDMDEDDDDHVDSDSIDSNDSDDDVDDTDDSHQ SDESHHSDESDELVTDFPTDLPATEVFTPVVPTVDTYDGRGDSVVYGL
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	SPP1
<b>Expression Region</b>	17-167aa
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	22.8 kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20 <sup>o</sup> /-80 <sup>o</sup> . Our default final concentration of glycerol is 50%. Customers could use it as reference.
-----------------------	---

<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.
-------------------	--



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.