




Recombinant Human Protein mono-ADP-ribosyltransferase PARP9 (PARP9)

Product Code	CSB-EP815575HU(A4)
Abbreviation	Recombinant Human PARP9 protein
Uniprot No.	Q8IXQ6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.103.
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MDFSMVAGAAAYNEKSGRITSLSLLFQKVFAQIFPQWRKGNTEECLPYKCSET GALGENYSWQIPINHNDFKILKNNERQLCEVLQNKFGCISTLVSPVQEGNSKSL QVFRKMLTPRIELSVWKDDLTHAVDAVVNAANEDLLHGGGLALALVKAGGFE IQEESKQFVARYGKVSAGEIAVTGAGRLPCKQIIHAVGPRWMEWDKQGCTGK LQRAIVSILNYVIYKNTHIKTVAIPALSSGIFQFPLNLCTKTIVETIRVSLQGKPM SNLKEIHLVSNEDPTVAAFKAASEFILGKSELGQETTPSFNAMVVNNTLQIVQ GHIEWQTADVIVNSVNPHDITVGPVAKSILQQAGVEMKSEFLATKAKQFQRSQ LVLVTKGFNLFCKYIYHVLWHSEFPKPKQILKHAMKECLEKCIQNITSISFPALGT GNMEIKKETAAEILFDEVLTFAKDHVKHQLTVKFVIFPTDLEIYKAFSSEMAKRS KMLSLNNYSVPQSTREEKRENGLEARSPAINLMGFNVEEMYEAHAWIQRILSL QNHIIENNHIYLGRKEHDILSQLQKTSSVSITEIISGRTELEIEGARADLIEVV MNIEDMLCKVQEEMARKKERGLWRSLGQWTIQQKTQDEMKENIIFLKCPVP PTQELLDQKKQFEKCGLQVLKVEKIDNEVLMAAFQRKKMMEEKLHRQPVSH RLFQQVPYQFCNVVCRVGFQRMYSTPCDPKYGAGIYFTKNLKNLAEKAKKISA ADKLIYVFEAEVLTGFFCQGHPLNIVPPPLSPGAIDGHDSVVDNVSSPETFVIFS GMQAIPQYLWTCTQEYVQSQDYSSGPMRPFQAHPWRGFASGSPVD
Research Area	Cancer
Source	E.coli
Target Names	PARP9
Expression Region	1-854aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	103.3 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.