




Recombinant Human Histone acetyltransferase type B catalytic subunit (HAT1)

Product Code	CSB-EP010143HU
Abbreviation	Recombinant Human HAT1 protein
Uniprot No.	O14929
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.153.
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AGFGAMEKFLVEYKSAVEKLAEYKCNTNTAIELKLVRFPELNDIRTFPEY THQLFGDDETAFGYKGLKILLYIAGSLSTMFRVEYASKVDENFDCVEADDVE GKIRQIIPPGFCTNTNDFLSLLEKEVDFKPFGLTLLHTYSVLSPTGGENFTFQIYK ADMTCRGFREYHERLQTFMLWFIETASFIDVDDERWHYFLVFEKYNKDGATLF ATVGYMTVYNYVYPDKTRPRVSQMLILTPFQGQGHGAQLLETVHRYYTEFP TVLDITAEDPSKSYVKLRDFVLVKLCQDLPCFSREKLMQGFNEDMVIEAQQKF KINKQHARRVYEILRLLVTDMSDAEQYRSYRLDIKRRLISPYKKKQRDLAKMRK CLRPEELTNQMNQIEISMQHEQLEESFQELVEDYRRVIERLAQE
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	HAT1
Expression Region	2-419aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	56.3 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.