




# Recombinant Human Elongation factor 1-alpha 1 (EEF1A1)(C31A)

<b>Product Code</b>	CSB-EP007409HU(M)
<b>Abbreviation</b>	Recombinant Human EEF1A1 protein (C31A)
<b>Uniprot No.</b>	P68104
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.159.
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MGKEKTHINIVVIGHVDSGKSTTTGHLIYKAGGIDKRTIEKFEKEAAEMGKGSFK YAWVLDKKAERERGITIDISLWKFFETSKYYYVTIIDAPGHRDFIKNMITGTSQAD CAVLIVAAGVGEFEAGISKNGQTREHALLAYTLGVKQLIVGVNKMDSTEPPYSQ KRYEEIVKEVSTYIKKIGYNPDTVAFVPISGWNGDNMLEPSANMPWFKGWKVT RKDGNASGTTLLEALDCILPPTTRPTDKPLRLPLQDVYKIGGIGTVPVGRVETGV LKPGMVVTFAPVNVTTTEVKSVEMHHEALSEALPGDNVGFNVKNVSVKDVRRG NVAGDSKNDPPMEAAGFTAQVIILNHPGQISAGYAPVLDCHTAHIACKFAELKE KIDRRSGKKLEDGPKFLKSGDAAIVDMVPGKPMCVESFSYDYPPLGRFAVRDM RQTVAVGVKAVDKKAAGAGKVTKSAQKAQKAK
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	E.coli
<b>Target Names</b>	EEF1A1
<b>Expression Region</b>	1-462aa(C31A)
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	63.1 kDa
<b>Protein Length</b>	Full Length

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.