




Recombinant Human Killer cell immunoglobulin-like receptor 3DL3 (KIR3DL3), partial, Biotinylated

Product Code	CSB-MP850818HUj8-B
Abbreviation	Recombinant Human KIR3DL3 protein, partial, Biotinylated
Uniprot No.	Q8N743
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.167.
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	QDKPFLSAWPGTVVSEGQHVTLQCRSRLGFNEFSLSKEDGMPVPELYNRIFR NSFLMGPVTPAHAGTYRCCSSHPHSPTGWSAPSNPVMVMTGVHRKPSLLAH PGPLVKSGETVILQCWSDVRFERFLLHREGITEDPLRLVGQLHDAGSQVNYSM GPMTPALAGTYRCFGSVTHLPYELSAPSDPLDIVVGLYGKPSLSAQPGPTVQ AGENVTLSCSSRSLFDIYHLSREAEAGELRLTAVLRVNGTFQANFPLGPVTHG GNYRCFGSFRALPHAWSDPDPLPVSVTGNSRHL
Research Area	Immunology
Source	Mammalian cell
Target Names	KIR3DL3
Expression Region	26-322aa
Tag Info	C-terminal hFc-Avi-tagged
Mol. Weight	62.9 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.