




Recombinant Human Tumor necrosis factor receptor superfamily member 10D (TNFRSF10D), partial

Product Code	CSB-EP887006HUa0
Abbreviation	Recombinant Human TNFRSF10D protein, partial
Uniprot No.	Q9UBN6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.178.
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	ATIPRQDEVPPQQTVPAPQQRRSLKEEECPAGSHRSEYTGACNPCTEGVDYTI ASNNLPSCLLCTVCKSGQTNKSSCTTTRDTCQCEKGSFQDKNSPEMCRCTCR TGCPRGMVKVSNCTPRSDIKCKNESAAASSTGKTPAAEETVTTLGMLASPYH
Research Area	Cancer
Source	E.coli
Target Names	TNFRSF10D
Expression Region	56-211aa
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.9 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.