



# Recombinant Human Flotillin-2 (FLOT2)

<b>Product Code</b>	CSB-EP619868HUa0
<b>Abbreviation</b>	Recombinant Human FLOT2 protein
<b>Uniprot No.</b>	Q14254
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.185.
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	GNCHTVGPNEALVVSGGCCGSDYKQYVFGGWAWAWWCISDTQRISLEIMTL QPRCEDVETAEGVALTVTGVAVQVKIMTEKELLAVACEQFLGKNVQDIKNVVLQ TLEGHLRSILGTLTVEQIYQDRDQFAKLVREVAAPDVGRMGIEILSFTIKDVYDK VDYLSSLGKTQTAVVQRDADIGVAEAERDAGIREAECKKEMLDVKFMADTKIA DSKRAFELQKSAFSEEVNIKTAEAQLAYELQGAREQQKIRQEEIEIEVVQRKKQI AVEAQEILRTDKELIATVRRPAEAEAHRIQQIAEGEKVKQVLLAQAEAEKIRKIG EAEAAVIEAMGKAEAEERMKLKAEAYQKYGDAKMALVLEALPQIAAKIAAPLTK VDEIVVLSGDNSKVTSEVNRLLAELPASVHALTGVDLSKIPLIKKATGVQV
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	FLOT2
<b>Expression Region</b>	2-428aa
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	52.9 kDa
<b>Protein Length</b>	Full Length of Mature Protein

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.