




Recombinant Mouse Kininogen-1 (Kng1)

Product Code	CSB-EP012479MO(F2)
Abbreviation	Recombinant Mouse Kng1 protein
Uniprot No.	O08677-2
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.186.
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	EEAQEIDCNDEAVFQAVDFSLKQFNPGVKSGNQYMLHRVIEGTKTDGSPTFY FKYLIKEGNCSAQSGLAWQDCDFKDAEEAATGECTATVGKREFFIVTQTCK IAPSKAPILKAYFPCIGCVHAISTDSPDLEPVLKHSIEHFNNTDHSFLTLRKVK SAHRQVVAGLNFDITYTIVQTNCSKERFPSLHGDCVALPNGDDGECRGNLFM DINNKIANFSQSCTLYSGDDLVEALPKPCPGCPRDIPVDSPELKEVLGHSIAQL NAENDHPFYKIDTVKKATSQVVAGTKYVIEFIARETKCSKESNTELAEDCEIKH LGQSLDCNANVYMRPWENKVVPTVKCQALDMTEMARRPPGFSPFRSVTVQE TKEGRTRLLRACEYKGRLSKAGAEPAPERQAESSQVKQ
Research Area	Signal Transduction
Source	E.coli
Target Names	Kng1
Expression Region	21-432aa
Tag Info	N-terminal 10xHis-tagged and C-terminal Flag-tagged
Mol. Weight	48.3 kDa
Protein Length	Full Length of Mature Protein of Isoform LMW

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.