




Recombinant Escherichia coli Glycine--tRNA ligase beta subunit (glyS)

Product Code	CSB-EP360602ENV
Abbreviation	Recombinant E.coli glyS protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P00961
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SEKTFLVEIGTEELPPKALRSLAESFAANFTAELDNAGLAHGTVQWFAAPRRLA LKVANLAEAQPDREIEIKRGPAIAQAFDAEGKPSKAAEGWARGCGITVDQAERL TTDKGEWLLYRAHVKGESTEALLPNMVATSLAKLPIPKLMRWGASDVHFVRPV HTVTLLLGDKVIPATILGIQSDRVIRGHRFMGEPEFTIDNADQYPEILRERGVIA DYEERKAKIKADAEAAARKIGGNADLSESLLEEVASLVEWPVVLTAKEEKFLA VPAEALVYTMKGDQKYFPVYANDGKLLPNFIFVANIESKDPQQIISGNEKVVRP RLADAEFFFNTDRKKRLEDNLPRLQTVLFQQQLGTLRDKTDRIQALAGWIAEQI GADVNHATRAGLLSKCDLMTNMVFETDTQGVMGMHYARHDGEAEDVAVAL NEQYQPRFAGDDLPSNPVACALAIADKMDTLGIFGIGQHPKGDKDPFALRRA ALGVLRIIVEKNLNDLQTLTEEAVRLYGDKLTNANVVDDVIDFMLGRFRAWYQ DEGYTVDTIQAVLARRPTRPADFDARMKAVSHFRTLDAALAAANKRVSNIL AKSDEVLSDRVNASTLKEPEEIKLAMQVVVLRDKLEPYFTEGGRYQDALVELAEL REPVDAFFDKVMVMVDDKELRINRLTMLEKLRELFRLVADISLLQ
Research Area	Others
Source	E.coli
Target Names	glyS
Expression Region	2-689aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	83.6 kDa

**Protein Length**

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.