




Recombinant Mouse Ras GTPase-activating protein-binding protein 1 (G3bp1)

Product Code	CSB-EP009116MO
Abbreviation	Recombinant Mouse G3bp1 protein
Uniprot No.	P97855
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVMEKPSPLLVGREFVRQYYTLLNQAPDMLHRFYGKNSSYAHGGLDSNGKP ADAVYGQKEIHRKVMSQNFTNCHTKIRHVDAHATLNDGVVVQVMGLLSNNNQ ALRRFMQTFVLAPEGSVANKFYVHNDIFRYQDEVFGGFVTEPQEESEEEVEEP EERQQTPEVVPDDSGTFYDQTVSNDLEEHLEEPVVEPEPEPEPEPEPEVSDI QEDKPEAALEEAAPDDVQKSTSPAPADVAPAQEDLRTFSWASVTSKNLPPSG AVPVTGTPPHVVKVPASQPRPESKPDSQIPPQRPQRDQRVREQRINIPPQRG PRPIREAGEPGDVEPRRMVRHPDSHQLFIGNLPHEVDKSELKDFQNFQGNVV ELRINSGGKLPNFGFVVFDDSEPVQKVLNRPIMFRGAVRLNVEEKKTRAARE GDRRDNRLRGPGGPRGGPSGGMRGPPRGGMVQKPGFGVGRGITTPRQ
Research Area	Cancer
Source	E.coli
Target Names	G3bp1
Expression Region	1-465aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	58.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.