




Recombinant Human Alpha-2-antiplasmin (SERPINF2)

Product Code	CSB-MP021085HU
Abbreviation	Recombinant Human SERPINF2 protein
Uniprot No.	P08697
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	MEPLGRQLTSGPNQEQVSPLTLLKLGNEPQQGQTALKSPPGVCSRDPDTPPEQT HRLARAMMAFTADLFLSLVAQTSTCPNLILSPLSVALALSHLALGAQNHTLQRLQ QVLHAGSGPCLPHLLSRLCQDLGPGAFRLAARMYLQKGFPIKEDFLEQSEQLF GAKPVSLTGKQEDDLANINQWVKEATEGKIQEFLSGLPEDTVLLLLLNAIHFQGF WRNKFDPSLTQRDSFHLDEQFTVPVEMMQARTYPLRWFLLEQPEIQVAHFPP KNNMSFVVLVPTHFEWNVSQVLANLSWDTLHPPLVWERPTKVRLPKLYLKHQ MDLVATLSQLGLQELFQAPDLRGISEQSLVSGVQHQSTLELSEVGVEAAAAT SIAMSRMSLSSFSVNRPFLLFFIFEDTTGLPLFVGSVRNPNPSAPRELKEQQDSP GNKDFLQSLKGFPRGDKLFGPDLKLVPPMEEDYPQFGSPK
Research Area	Cardiovascular
Source	Mammalian cell
Target Names	SERPINF2
Expression Region	28-491aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	53.1 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.