




Recombinant Human 2-5A-dependent ribonuclease (RNASEL), partial

Product Code	CSB-BP313223HU2
Abbreviation	Recombinant Human RNASEL protein, partial
Uniprot No.	Q05823
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	AAVEDNHLLIKAVQNEVDLVQQLLEGGANVNFQEEEGGWTPHNAVQMSRE DIVELLRRHGADPVLRRKNGATPFILAAIAGSVKLLKFLSKGADVNECDFYGF AFMEAAVYGKVKALKFLYKRGANVNLRRKTKEDQERLRKGGATALMDAAEK HVEVLKILLDEMADVACDNMGRNALIHALLSSDDSDVEAITHLLLDHGADV VRGERGKTPLILAVEKKHLGLVQRLLLEQEHIEINDTSDGKTALLLAVELKLLKIA ELLCKRGASTDCGDLVMTARRNYDHSLVKVLVSHGAKEDFHPPAEDWKPQSS HWGAALKDLHRIYRPMIGKLFIDEKYKIADTSEGGIYLGfyEKQEVAVKTFCE GSPRAQREVSLQSSRENSHLVTFYGSSEHRGHLFVCVTLCEQTLEACLDVH RGEDVENEDEFARNVLSSIFKAVQELHLSCGYTHQDLQPQNILIDSKAAHLA DFDKSIKWAGDPQEVKRDLEDLGRLLVLYVVKKGSISFEDLKAQSNEEVQLSP DEETKDLIHRFHPGEHVRDCLSDLLGHPFFWTWESRYRTLNRVGNESDIKTR KSESEILRLLQPGPSEHSKSFDKWTTKINECVMMKMNKfyEKRGNFYQNTVG DLLKFIRNLGEHIDEKHKMKLKGDPsLYFQKTFPDLVIYVYTKLQNTeyRKH FPQTHSPNK
Research Area	Epigenetics and Nuclear Signaling
Source	Baculovirus
Target Names	RNASEL
Expression Region	21-724aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	81.3 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.