



AMHR2 (Murlentamab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA614985MB1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16671
Immunogen	Recombinant Human AMHR2 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	3-4 weeks
Target Names	AMHR2
Notes	<p>Validation Status</p> <p>Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.</p> <p>Guaranteed Quality</p> <p>① Antibody purity > 95% tested by SDS-PAGE.</p> <p>② Endotoxin level < 0.1EU/ug tested by LAL method.</p>