




Recombinant Human AP-1 complex subunit beta-1 (AP1B1), partial

Product Code	CSB-EP604616HUc7
Abbreviation	Recombinant Human AP1B1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	Q10567
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MTDSKYFTTTKKGEIFELKAELNSDKKEKKKEAVKKVIASMTVGKDVSAALFPDV VNCMQTDNLELKKLVYLYLMNYAKSQPDMAIMAVNTFVKDCEDPNPLIRALAV RTMGCIRVDKITEYLCEPLRKCLKDEDPYVRKTAAVCVAKLHDINAQLVEDQGF LDTLKDLSNSPMVVANAVAALSEIAESHPSNLLDLNPQSINKLLTALNECTE WGQIFILDCLANYMPKDDREAQSICERVTPRLSHANSVVLSAVKVLKFMEM LSKDLDYGGTLLKKLAPPLVTLSSAEPQLQYVALRNINLIVQKRPEILKHEMKVFF VKYNDPIYVKLEKLDIMIRLASQANIAQVLAELKEYATEVDVDFVRKAVRAIGRC AIKVEQSAERCVSTLLDLIQTKVNYVVQEAIVVIKDFRKYPNKYESVIATLCENL DSLDEPEARAAAMIWIVGEYAERIDNADELLESFLEGFHDESTQVQLQLLTAIVKL FLKKPTETQELVQQVLSLATQSDNPDLRDRGYIYWRLSTDPVAAKEVVLAE KPLISEETDLIEPTLLDELICYIGTLASVYHKPPSAFVEG
Research Area	Signal Transduction
Source	E.coli
Target Names	AP1B1
Expression Region	1-584aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	72.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.