



Recombinant Human Cadherin-3 (CDH3), partial

Product Code	CSB-MP005052HU2
Abbreviation	Recombinant Human CDH3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	P22223
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	EPCRAVFREA EVTLEAGGAEQEPGQALGKVFMGCPGQEPALFSTDNDDFTV RNGETVQERRSLKERNPLKIFPSKRILRRHKRDWV VAPISVPENKGPFPQRL NQLKSNKDRDTKIFY SITGPGADSPPEGVFAVEKETGWLLLNKPLDREEIAKYE LFGHAVSENGASVEDPMNISIIVTDQNDHKPKFTQDTRGVS VLEGVLPGT SVM QVTATDEDDAIYTYNGVVAYSIH SQEPKDPHDL MFTIHRSTGTISVISSGLDREK VPEYTLTIQATDMDGDGSTTTAVAVVEILDANDNAPMFD P QKYEAHV PENAVG HEVQRLTVTDLDAPNSPAWRATYLIMGGDDGDHFTITTHPESNQGILTTRKGL DFEAKNQHTLYVEVTNEAPFVLKLP TSTATIVVHVEDVNEAPV FVPPSKVVEVQ EGIPTGEPVCVYTAEDPDKENQKISYRILRDPAGWLAMD PDSGQVTAVGTLDR EDEQFVRNNIYEVMLAMDNGSPPTTGTG TLLL TLIDVNDHGPVPEPRQITICN QSPVRQVLNITDKDLS PHTSPFQAQLTDDSDIYWTAEVNEEGDTVVLSLKKFLK QDTYDVHLSLSDHGNKEQLTVIRATVCDCHGHVETCPGPWKGG
Research Area	Cancer
Source	Mammalian cell
Target Names	CDH3
Expression Region	25-654aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	70.8 kDa
Protein Length	Partial

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.