



Recombinant Mouse Vascular endothelial growth factor receptor 2 (Kdr), partial

Product Code	CSB-MP012145MO
Abbreviation	Recombinant Mouse Kdr protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	P35918
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	ASVGLPGDFLHPPKLSTQKDILTILANTTLQITCRGQRDLDWLWPNAQRDSEE RVLVTECGGGDSIFCKTLTIPRVVGNDTGAYKCSYRDVDIASTVYVYVRDYRS PFIASVSDQHGIVYITENKNKTVVPCRGSSISNLNVSLCARYPEKRFVDPGNRIS WDSEIGFTLPSYMISYAGMVFCEAKINDETYQSIMYIVVVVGYRIYDVLSPPEI ELSAGEKLVLNCTARTELVGLDFTWHSPPSKSHHKIVNRDVKPFPGTVAKM FLSTLTIESVTKSDQGEYTCVASSGRMIKRNRTFVRVHTKPFIAFGSGMKSLVE ATVGSQVRIPVKYLSYPAPDIKWYRNGRPIESNYTMIVGDELTIMEVTERDAGN YTVILTNPISMEKQSHMVSLVNVNPPQIGEKALISPMSYQYGTMTLTCTVYA NPPLHHIQWYWQLEEACSYRPGQTSPYACKEWRHVEDFQGGNKIEVTKNQY ALIEGKNKTVSTLVIQAANVSALYKCEAINKAGRGERSVIFHVIRGPEITVQPAA QPTEQESVSLCTADRNTFFENLTWYKLGSAQTSVHMGESLTPVCKNLDALWK LNGTMFSNSTNDILIVAFQNASLQDQGDYVCSAQDKKTKKRHCLVKQLIILERM APMITGNLENQTTTIGETIEVTCPASGNPTPHITWFKDNETLVEDSGIVLRDGN RNLTIRVRKEDGGLYTCQACNVLGCARAETLFIIEGAQEKTNLE
Research Area	Others
Source	Mammalian cell
Target Names	Kdr
Expression Region	20-762aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	111.9 kDa

**Protein Length**

Partial

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.