




# Recombinant Human Adenosine receptor A2a (ADORA2A) (F168A)-Detergent

<b>Product Code</b>	CSB-MP001376HU(A4)(M)-D
<b>Abbreviation</b>	Recombinant Human ADORA2A protein (F168A)-Detergent
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
<b>Uniprot No.</b>	P29274
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Protein
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MPIMGSSVYITVELAIAVLAILGNVLCWAVWLNSNLQNVTNYFVVSLLAAADIAV GVLAIPTAITISTGFCAACHGCLFIACFVLVLTQSSIFSLLAIAIDRYIAIRIPLRYNG LVTGTRAKGIIAICWVLSFAIGLTPMLGWNNCGQPKEGKNHSQGCQEGQVACL AEDVVPMNYMVYFNFFACVLVPLLLMLGVYLRIFLAARRQLKQMESQPLPGER ARSTLQKEVHAAKSLAIIIVGLFALCWLP LHIINCFTFFCPDCSHAPLWLMYLAIVL SHTNSVVPFIYAYRIREFRQTRKIIIRSHVLRQQEPFKAAGTSARVLAAHGSD GEQVSLRLNGHPPGVWANGSAPHPERRPNGYALGLVSGGSAQESQGNTGL PDVELLSHELKGVCPPEPGLDDPLAQQDGAGVS
<b>Research Area</b>	Neuroscience
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ADORA2A
<b>Expression Region</b>	1-412aa(F168A)
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-Flag-tagged and C-terminal Twin-Strep-tagged
<b>Mol. Weight</b>	50.2 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	 <p>(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.</p>



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.