



# Recombinant Mouse Cannabinoid receptor 1 (Cnr1)-VLPs

<b>Product Code</b>	CSB-MP005678MO
<b>Abbreviation</b>	Recombinant Mouse Cnr1 protein-VLPs
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
<b>Uniprot No.</b>	P47746
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
<b>Product Type</b>	Recombinant Protein
<b>Sensitivity</b>	Not Test
<b>Purity</b>	The purity information is not available for VLPs proteins.
<b>Sequence</b>	MKSILDGLADTTFRITITDALLYVGSNDIQYEDIKGDMSKLGYPQKFPLTSFRG SPFQEKMTAGDNSPLVPAGDTTINITEFYNKSLSSFKENEDNIQCGENFMDME CFMILNPSQQLAIAVLSLTLGTFTVLENLLVLCVILHSRSLRCRPSYHFIGSLAVA DLLGSVIVVYSFVDFHVFHRKDSPNVFLFKLGGVTASFTASVGSFLTAIDRYISI HRPLAYKRIVTRPKAVVAFCLMWTAIAVIAVLP LLGW NCKKLQSVCSDFIPLIDET YLMFWIGVTSVLLLFIVYAYMYILWKAHSHAVRMIQRGTQKSIIIHTSEDGKVQV TRPDQARMDIRLAKTLVLILVLIICWGPLLAIMVYDVF GKMNKLIKTVFAFCSML CLLNSTVNPIIYALRSKDLRHAFRSMFPSCGTAQPLDNMGS DCLHKHANN TASMHRAAESCICKSTVKIAKVTMSVSTDTSAEAL
<b>Research Area</b>	Cancer
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Cnr1
<b>Expression Region</b>	1-473aa
<b>Notes</b>	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20?/-80? upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
<b>Tag Info</b>	C-terminal 10xHis-tagged
<b>Mol. Weight</b>	54.3 kDa
<b>Protein Length</b>	Full Length

**Image**


CSB-MP005678MO is detected by Mouse



anti-6\*His monoclonal antibody

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.