




Recombinant Human Tribbles homolog 3 (TRIB3)

Product Code	CSB-EP847220HUe1
Abbreviation	Recombinant Human TRIB3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	Q96RU7
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	MRATPLAAPAGSLSRKKRLELDDNLDTERPVQKRARSGPQPRLPPCLLPLSPP TAPDRATAVATASRLGPYVLEPEEGGRAYQALHCPTGTEYTCKVYPVQEALA VLEPYARLPPHKHVARPTEVLAGTQLLYAFFTRTHGDMHSLVRSRHRPEPEA AVLFRQMATALAHCHQHGLVLRDLKLCRFVFADRERKKLVLENLEDSCVLTGP DDSLWDKHACPAYVGPEILSSRASYSKGAADVWSLGVLFMLAGHYPFQDS EPVLLFGKIRRGAYALPAGLSAPARCLVRCLLRREPAERLTATGILLHPWLRQD PMPLAPTRSHLWEAAQVVPDGLGLDEAREEEGDREVVLYG
Research Area	Cancer
Source	E.coli
Target Names	TRIB3
Expression Region	1-358aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	Tag-Free
Mol. Weight	39.7 kDa
Protein Length	Full Length
Image	 <p>(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.</p>
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°/-80°. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.