




Recombinant Human Delta-type opioid receptor (OPRD1)

Product Code	CSB-CF016358HUm12
Abbreviation	Recombinant Human OPRD1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	P41143
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Sensitivity	Not Test
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MEPAPSAGAE LQPPLFANASDAYPSACPSAGANASGPPGARSASSLALAIIT ALYSAVCAVGLLG NVLVMFGIVRYTKMKTATNIYIFNLALADALATSTLPFQSAK YLMETWPFGE LLCKAVLSIDYYNMF TSIFTLTMM SVDRYIAVCHPVKALDFRTP AKAKLINICIWVLASGVGVPIMVMAVTRPRDGAVVCMLQFPSPSWYWDTVTKI CVFLFAFVVPIIITVCYGLMLLR LRSVRLLSGSKEKDRSLRRITRMVLVVVGAFV VCWAPIHIFVIVWTLVDIDRRDPLVVAALHLCIALGYANSSLNPVLYAFLDENFK RCFRQLCRKPCGRPD PSSF SRAREATARERV TACTPSDGPGGGAAA
Research Area	Others
Source	in vitro E.coli expression system
Target Names	OPRD1
Expression Region	1-372aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Twin-Strep-tagged
Mol. Weight	44.9 kDa
Protein Length	Full Length
Image	 <p>(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.</p>
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.