




HOXB9 Recombinant Monoclonal Antibody

Product Code	CSB-RA184200A0HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17482
Immunogen	A synthesized peptide derived from Human HOXB9 protein
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HOXB9
Clone No.	14G3

Image



Western Blot

Positive WB detected in: A549 whole cell lysate(30µg), Jurkat whole cell lysate(30µg), U251 whole cell lysate(30µg), COLO205 whole cell lysate(30µg)
 All lanes: HOXB9 antibody at 1:1000
 Secondary
 Goat polyclonal to rabbit IgG at 1/40000 dilution
 Predicted band size: 28 kDa
 Observed band size: 28 kDa
 Exposure time?2min

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.