



EPCAM (Tucotuzumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA007717MB3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16422
Immunogen	Recombinant Human EPCAM protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG2 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	EPCAM
Accession NO.	1E1

Image


The Binding Activity of Human EPCAM with Anti-EPCAM recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human EPCAM (CSB-MP007717HU-B) at 2 µg/mL can bind Anti-IFNG recombinant antibody. The EC₅₀ is 1.357 - 1.493 ng/mL.