




FCRL5 Recombinant Monoclonal Antibody


Product Code	CSB-RA853494MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD9
Immunogen	Recombinant Human FCRL5 protein
Species Reactivity	Human
Tested Applications	ELISA, SPR; Recommended dilution: ELISA:1:1000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCRL5
Clone No.	17G2

Image



The Binding Activity of Human FCRL5 with Anti-FCRL5 Recombinant Antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Human FCRL5 (CSB-MP853494HU1) at 2 µg/mL can bind Anti-FCRL5 recombinant antibody. The EC50 is 1.107-1.253 ng/mL.



The Binding Activity of Human FCRL5 with Anti-FCRL5 Recombinant Antibody

Activity: FCRL5 Recombinant Monoclonal Antibody captured on Protein A Chip can bind Recombinant Human FCRL5 (CSB-MP853494HU1) with an affinity constant of 0.436 nM as detected by MetaSPR Assay (WeSPRTM 200).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.