




CDH3 Recombinant Monoclonal Antibody

Product Code	CSB-RA005052MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22223
Immunogen	Recombinant Human CDH3 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:1000-1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CDH3
Clone No.	1F9

Image



The Binding Activity of Human CDH3 with Anti-CDH3 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CDH3 (CSB-MP005052HU7) at 2 µg/mL can bind Anti-CDH3 recombinant antibody. The EC₅₀ is 2.147-2.507 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.