





# SNCA Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA021912MA2HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P37840
<b>Immunogen</b>	Recombinant Human SNCA protein
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, SPR; Recommended dilution: ELISA:1:1000-1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	hIgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Recombinant Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	SNCA
<b>Clone No.</b>	20F1

## Image



The Binding Activity of Human SNCA with Anti-SNCA recombinant antibody  
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SNCA(CSB-MP021912HU) at 2 µg/mL can bind anti-SNCA recombinant antibody. The EC50 is 3.249-7.254 ng/mL.



The Binding Activity of Human SNCA with Anti-SNCA recombinant antibody  
Activity: SNCA Recombinant Monoclonal Antibody captured on Protein A Chip can bind Human SNCA(CSB-MP021912HU) with an affinity constant of 19.7 nM as detected by MetaSPR Assay (WeSPR™ 200).

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.