




EGFR Recombinant Monoclonal Antibody

| | |
|----------------------------|--|
| Product Code | CSB-RA159341MA2HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P00533 |
| Immunogen | Recombinant Human EGFR protein |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, FC; Recommended dilution: FC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Affinity-chromatography |
| Isotype | mIgG2a |
| Clonality | Monoclonal |
| Product Type | Recombinant Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | EGFR |
| Clone No. | 10F6 |

Image



Untransfected CHO-K1*10⁶ cells), washed and then followed by APC-conjugated anti-Mouse IgG Fc antibody and analyzed with flow cytometry.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.