



CD274 Recombinant Monoclonal Antibody

Product Code	CSB-RA878942MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZQ7
Species Reactivity	Mouse
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG4(S228P) kappa
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	Cd274
Clone No.	368A.4H1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.