



Ctla4 Recombinant Monoclonal Antibody

Product Code	CSB-RA006163MA1MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09793
Species Reactivity	Mouse
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	scFv-kappa-heavy
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Mus musculus (Mouse)
Target Names	CTLA4
Clone No.	9D9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.