



LRRC32 (Livmoniplimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA618886MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14392
Immunogen	Recombinant Human LRRC32 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG4(S228P) kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRRC32
Accession NO.	7A9