





CD200 (Samalizumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA004895MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41217
Immunogen	Recombinant Human CD200 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD200
Accession NO.	7B10

Image



The Binding Activity of Human CD200 with Anti-CD200 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD200 (CSB-MP004895HU1) at 2 µg/mL can bind Anti-CD200 recombinant antibody. The EC₅₀ is 2.984 - 3.700 ng/mL.



The purity of CD200 was greater than 95% as determined by SEC-HPLC