



# IL-31 (Lokivetmab Biosimilar) Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA984018MA1DO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	C7G0W1
<b>Immunogen</b>	Recombinant Human IL31 (Canine) protein
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, FACS
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	Human IgG1 kappa
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Biosimilar Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	IL31 (Canine)
<b>Accession NO.</b>	7F2