





SELPLG (Neihulizumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA621767MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14242
Immunogen	Recombinant Human CD162 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 lambda2
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD162
Accession NO.	7H4

Image



The Binding Activity of Human SELPLG with Anti-SELPLG recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SELPLG (CSB-MP621767HU1) at 2 µg/mL can bind Anti-SELPLG recombinant antibody. The EC₅₀ is 2.709 - 3.123 ng/mL.



The purity of CD162 was greater than 95% as determined by SEC-HPLC