





IL4R (Manfidokimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA011661MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24394
Immunogen	Recombinant Human IL4R protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL4R
Accession NO.	6B2

Image



The Binding Activity of Human IL4R with Anti-IL4R recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IL4R (CSB-MP011661HU1) at 2 µg/mL can bind Anti-IL4R recombinant antibody. The EC₅₀ is 2.712 - 2.990 ng/mL.



The purity of IL4R was greater than 95% as determined by SEC-HPLC