





TGM2 (Zampilimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA023462MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21980
Immunogen	Recombinant Human TGM2 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 lambda1
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TGM2
Accession NO.	6D4

Image



The Binding Activity of Human TGM2 with Anti-TGM2 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TGM2 (CSB-EP023462HU) at 2 µg/mL can bind Anti-TGM2 recombinant antibody. The EC₅₀ is 3.135 - 3.679 ng/mL.



The purity of TGM2 was greater than 95% as determined by SEC-HPLC