





IL21 (Avizakimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA872537MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBE4
Immunogen	Recombinant Human IL21 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL21
Accession NO.	6E5

Image



The Binding Activity of Human IL21 with Anti-IL21 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IL21 (CSB-EP872537HU) at 2 µg/mL can bind Anti-IL21 recombinant antibody. The EC₅₀ is 24.61 - 29.23 ng/mL.



The purity of IL21 was greater than 95% as determined by SEC-HPLC