





CD22 (Moxetumomab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA004900MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20273
Immunogen	Recombinant Human CD22 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD22
Accession NO.	6G7

Image



The Binding Activity of Human CD22 with Anti-CD22 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD22 (CSB-MP004900HU) at 1 µg/mL can bind Anti-CD22 recombinant antibody. The EC₅₀ is 0.5197-0.6255 ng/mL.



The purity of CD22 was greater than 95% as determined by SEC-HPLC