





IFNG (Fontolizumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA011050MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01579
Immunogen	Recombinant Human IFNG protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IFNG
Accession NO.	1A9

Image



The Binding Activity of Human IFNG with Anti-IFNG recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IFNG (CSB-YP011050HU1) at 2 µg/mL can bind Anti-IFNG recombinant antibody. The EC₅₀ is 1.123 - 1.230 ng/mL.



The purity of IFNG was greater than 95% as determined by SEC-HPLC