





CD52 (Gatralimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA004943MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31358
Immunogen	Recombinant Human CD52 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG2 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD52
Accession NO.	1F2

Image



The Binding Activity of Human CD52 with Anti-CD52 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD52 (CSB-MP004943HUd7) at 2 µg/mL can bind Anti-CD52 recombinant antibody. The EC₅₀ is 9.793 - 12.60 ng/mL.



The purity of CD52 was greater than 95% as determined by SEC-HPLC