



Toxin A (Tosatoxumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA357716MA1FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09616
Immunogen	Recombinant Staphylococcus aureus Toxin A protein
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Staphylococcus aureus
Target Names	Staphylococcus aureus Toxin A
Accession NO.	4C11