



# Toxin A (Tosatoxumab Biosimilar) Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA357716MB1FKZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P09616
<b>Immunogen</b>	Recombinant Staphylococcus aureus Toxin A protein
<b>Species Reactivity</b>	Staphylococcus aureus
<b>Tested Applications</b>	ELISA, FACS
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Isotype</b>	Human IgG1 kappa
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Biosimilar Antibody
<b>Immunogen Species</b>	Staphylococcus aureus
<b>Target Names</b>	Staphylococcus aureus Toxin A
<b>Accession NO.</b>	4C11