





LIF (Falbikitug Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA012928MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15018
Immunogen	Recombinant Human LIF protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LIF
Accession NO.	4D12

Image



The Binding Activity of Human LIF with Anti-LIF recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human LIF (CSB-MP012928HU) at 2 µg/mL can bind Anti-LIF recombinant antibody. The EC₅₀ is 9.793 - 12.60 ng/mL.



The purity of LIF was greater than 95% as determined by SEC-HPLC