



SLITRK6 (Sirtratumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA867155MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5Y7
Immunogen	Recombinant Human SLITRK6 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLITRK6
Accession NO.	4G3