





MUC16 (Sofituzumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA704410MA3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI7
Immunogen	Recombinant Human MUC16 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG3 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MUC16
Accession NO.	4H4

Image



The Binding Activity of Human MUC16 with Anti-MUC16 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human MUC16 (CSB-MP704410HU3c7) at 2 µg/mL can bind Anti-MUC16 recombinant antibody. The EC₅₀ is 8.605 - 10.18ng/mL.



The purity of MUC16 was greater than 95% as determined by SEC-HPLC