



VTCN1 (Puxitatug Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA801822MB3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z7D3
Immunogen	Recombinant Human VTCN1 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VTCN1
Accession NO.	5A5

Image



The Binding Activity of Human VTCN1 with Anti-VTCN1 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human VTCN1 (CSB-MP801822HU1) at 2 µg/mL can bind Anti-VTCN1 recombinant antibody. The EC₅₀ is 8.605 - 10.18ng/mL.



The purity of VTCN1 was greater than 95% as determined by SEC-HPLC