



FCRL5 (DFRF4539A (RG-7598) Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA853494MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD9
Immunogen	Recombinant Human FCRL5 protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCRL5
Accession NO.	3F6