



# GD3 ganglioside (Ecromeximab Biosimilar) Recombinant Monoclonal Antibody

<b>Product Code</b>	CSB-RA044639I1m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	GD3 ganglioside protein
<b>Species Reactivity</b>	All
<b>Tested Applications</b>	ELISA, FACS
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	0.01M PBS, pH 7.4
<b>Purification Method</b>	Affinity-chromatography
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Biosimilar Antibody
<b>Immunogen Species</b>	All
<b>Target Names</b>	GD3 ganglioside
<b>Accession NO.</b>	4A9