



GD3 ganglioside (Ecromeximab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA044639MB1I
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	GD3 ganglioside protein
Species Reactivity	All
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	All
Target Names	GD3 ganglioside
Accession NO.	4A9