



# Recombinant Clostridium perfringens Perfringolysin O (pfo)

<b>Product Code</b>	CSB-EP314820CMBb0
<b>Abbreviation</b>	Recombinant Clostridium perfringens pfo protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
<b>Uniprot No.</b>	P0C2E9
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium perfringens (strain 13 / Type A)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	≥95% as determined by SDS-PAGE.
<b>Sequence</b>	KDITDKNQSIDSGISSLSYNRNEVLASNGDKIESFVPKEGKKTGNKFIVVERQK RSLTTSPVDISIIDS VNDRTYPGALQLADKAFVENRPTILMVKRKPININIDLPLGLK GENSIKVDDPTYGKVS GAIDELVSKWNEKYSSTHTLPARTQYSESMVYSKSQI SSALNVNAKVLENSLGVDFNAVANNEKKVMILAYKQIFYTVSADLPKNPSDLFD DSVTFNDLKQKGVSN EAPPLMVS NVAYGRTIYVKLETTSSSKDVQAAFKALIKN TDIKNSQQYKDIYENSSFTAVVLGGDAQEHNKVVTKDFDEIRKVIKDNATFSTK NPAYPISYTSVFLKDNSVA AVHNKTDYIETTSTEYSK GKINLDHSGAYVAQFEV AWDEVSYDKEGNEVLTHKTWDGNYQDKTAHYSTVIPLEANARNIRIKARECTG LAWEWWRD VISEYDVPLTNNIN VSIWGTTL YPGSSITYN
<b>Research Area</b>	Microbiology
<b>Source</b>	E.coli
<b>Target Names</b>	pfo
<b>Expression Region</b>	29-500aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	56.3 kDa
<b>Protein Length</b>	Full Length of Mature Protein

**Image**


(Tris-Glycine gel) Discontinuous SDS-PAGE  
(reduced) with 5% enrichment gel and 15%



separation gel.

---

**Typical**

No

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.