




# Recombinant Human Beta-klotho (KLB), partial

<b>Product Code</b>	CSB-MP769815HU
<b>Abbreviation</b>	Recombinant Human KLB protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
<b>Uniprot No.</b>	Q86Z14
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	≥95% as determined by SDS-PAGE.
<b>Sequence</b>	FSGDGRAIWSKNPNFTPVNESQLFLYDTFPKNFFWGIGTGALQVEGSWKKDG KGPSIWDHFIHTHLKNVSSTNGSSDSYIFLEKDLSALDFIGVSFYQFSISWPRLF PDGIVTVANAKGLQYYSTLLDALVLRNIEPIVTLYHWDLPLALQEKEYGGWKNNDT IIDIFNDYATYCFQMFGDRVKYWITIHPYLVAWHGYGTGMHAPGEKGNLAAV YTVGHNLKAHASKVWHNYNTHFRPHQKGWLSITLGSWHIEPNRSENTMDIFKC QQSMVSVLGWFANPIHGDGDYPEGMRKKLFSVLPFSEAEKHEMRGTADFFA FSFGPNNFKPLNTMAKMGQNVSLNLREALNWIKLEYNNPRILIAENGWFTDSR VKTEDTTAIYMMKNFLSQVLQAIRLDEIRVFGYTAWSLLDGFEWQDAYTIRRGL FYVDFNSKQKERKPKSSAHYYKQIIRENGFSLKESTPDVQGGQFPDFSWGVT ESVLKPESVASSPQFSDPHLYVWNATGNRLLHRVEGVRLKTRPAQCTDFVNIK KQLEMLARMKVTHYRFALDWASVLPTGNLSAVNRQALRYRCVVSEGLKLG SAMVTLYYPHTAHLGLPEPLLHADGWLNPSTAEAFQAYAGLCFQELGDLVKL WITINEPNRLSDIYNRSGNDTYGAAHLLVAHALAWRLYDRQFRPSQRGAVSL SLHADWAEPANPYADSHWRAAERFLQFEIAWFAEPLFKTGDYPAAMREYIAS KHRRGLSSSALPRLTEAERRLLKGTVDFCALNHFTTRFVMHEQLAGSRYDSD RDIQFLQDITRLSSPTRLAVIPWGVRRKLLRWVRRNYGDMDIYITASGIDDQALE DDRLRKYLLGKYLQEVKAYLIDKVRKIGYYAFKLAEEKSKPRFGFFTSDFKAK SSIQFYNKVISSRGFPFENSSSRCSQTQENT
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KLB
<b>Expression Region</b>	53-983aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.



<b>Tag Info</b>	C-terminal 10xHis-tagged
<b>Mol. Weight</b>	108.8 kDa
<b>Protein Length</b>	Partial
<b>Image</b>	 (Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.
<b>Typical</b>	No
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.